

**AERB/RSD/PET/QAR/Rev.3/2021**

**QUALITY ASSURANCE TESTS FOR  
PET-CT / PET-MR/ PET IMAGING EQUIPMENT**

**RADIOLOGICAL SAFETY DIVISION  
ATOMIC ENERGY REGULATORY BOARD,  
NIYAMAK BHAVAN, ANUSHAKTINAGAR  
MUMBAI-400 094**

## Quality Assurance Tests for PET-CT/ PET-MR/PET Imaging Equipment

Institute Name	
Institute number in eLORA	
Installation ID in eLORA	
<b>Details of PET-CT / PET-MR/ PET</b>	
Type of Equipment(PET-CT / PET-MR/ PET)	
Manufacturer of the equipment	
Model name of the equipment	
The type/name of crystal used	
The number of crystals	
The number of rings of detectors	
<b>CT component details of PET-CT *</b>	
Manufacturer of the CT component	
Model name of the CT component	
Type approval number of the CT component	

\*Mention NA if CT component is not available

### I. Test related to PET component

#### 1. Spatial resolution (FWHM & FWTM) in 2D and 3D modes (as applicable\*) at the centre:-

Description	2 D Mode			3D Mode		
	Measured (mm)	Reference (mm)	Tolerance	Measured (mm)	Reference (mm)	Tolerance
FWHM in Tangential Direction			≤ 1.05 of Reference			≤ 1.05 of Reference
FWHM in Radial Direction			≤ 1.05 of Reference			≤ 1.05 of Reference
FWHM in Axial Direction			≤ 1.05 of Reference			≤ 1.05 of Reference
FWTM in Tangential Direction			≤ 1.05 of Reference			≤ 1.05 of Reference
FWTM in Radial Direction			≤ 1.05 of Reference			≤ 1.05 of Reference
FWTM in Axial Direction			≤ 1.05 of Reference			≤ 1.05 of Reference

\* Please indicate NA in case a particular mode is not applicable

Note: The ratio between FWTM and FWHM should be in the range 1.8–2.0 if the manufacturer reference value is not available for the same

**2. Spatial resolution (FWHM & FWTM) in 2D and 3D modes (as applicable\*) at 10cm from the central axis:-**

Description	2 D Mode			3D Mode		
	Measured (mm)	Reference (mm)	Tolerance	Measured (mm)	Reference (mm)	Tolerance
FWHM in Tangential Direction			≤ 1.05 of Reference			≤ 1.05 of Reference
FWHM in Radial Direction			≤ 1.05 of Reference			≤ 1.05 of Reference
FWHM in Axial Direction			≤ 1.05 of Reference			≤ 1.05 of Reference
FWTM in Tangential Direction			≤ 1.05 of Reference			≤ 1.05 of Reference
FWTM in Radial Direction			≤ 1.05 of Reference			≤ 1.05 of Reference
FWTM in Axial Direction			≤ 1.05 of Reference			≤ 1.05 of Reference

\* Please indicate NA in case a particular mode is not applicable

Note: The ratio between FWTM and FWHM should be in the range 1.8–2.0 if the manufacturer reference value is not available for the same

**3. Sensitivity in 2D and 3D modes (as applicable\*):-**

Description	Measured	Reference	Tolerance
<b>Sensitivity</b> in CPS/(Bq/ml) for 2 D Mode			≥ 0.95 of Reference
<b>Sensitivity</b> in CPS/(Bq/ml) for 3 D Mode			≥ 0.95 of Reference

\* Please indicate NA in case a particular mode is not applicable

**4. Scatter Fraction in 2D and 3D modes (as applicable\*):-**

Description	Measured	Reference	Tolerance
Scatter Fraction (2 D Mode) in %			≤ 1.05 of Reference
Scatter Fraction (3 D Mode) in %			≤ 1.05 of Reference

\* Please indicate NA in case a particular mode is not applicable

**5. Noise Equivalent Count Rate (NECR) in 2D and 3D modes (as applicable\*):-**

Description	2 D Mode			3D Mode		
	Measured	Reference	Tolerance	Measured	Reference	Tolerance
Peak NECR (kcps) in 2 D Mode			$\geq 0.95$ of Reference			$\geq 0.95$ of Reference
Activity (MBq) at which Peak NECR is achieved in 2D mode			$\geq 0.95$ of Reference			$\geq 0.95$ of Reference
Peak NECR (kcps) in 3 D Mode			$\geq 0.95$ of Reference			$\geq 0.95$ of Reference
Activity (MBq) at which Peak NECR is achieved in 3D mode			$\geq 0.95$ of Reference			$\geq 0.95$ of Reference

**6. Energy resolution:-**

Description	Measured	Reference	Tolerance
Energy resolution			$\leq 1.05$ of Reference

**7. Image quality, accuracy of attenuation & scatter correction and accuracy of radioactivity quantitation:-**

Description	Measured	Reference	Tolerance
<b>Image quality</b>			
The percent contrast for each hot sphere			$\pm 5\%$ of Reference
The percent contrast for each cold sphere			$\pm 5\%$ of Reference
The percent background variability			$\pm 5\%$ of Reference
<b>Accuracy of attenuation and scatter corrections</b>			
Percent relative error for each slice in lung slice			$\pm 5\%$ of Reference
<b>Accuracy of radioactivity quantitation</b>			
Radioactivity concentration in MBq/ml			$\pm 5\%$ of Reference

**8. Coincidence timing resolution for TOF positron emission tomography (as applicable\*):-**

Description	Measured	Reference	Tolerance
Timing FWHM (Resolution Timing)			≤ 1.05 of Reference

\* Please indicate NA in case not applicable

**9. Non-Uniformity of the reconstructed image:-**

Description	Measured	Reference	Tolerance
Mean non-uniformity (% NU)			≤ 1.05 of Reference
Coefficient of variation			
Standard Deviation (SD)			

**10. Normalization: -**

Description	Observation	Tolerance
Whether uniformity observed in the normalization sinogram	Yes/No	Observance of Uniformity in the sinogram

**11. Accuracy of PET/CT image registration:-**

Description	Measured	Reference	Tolerance
Registration Value for 512x512 matrix			±1 pixel (or ±1 mm, whichever is smaller)

**UNDERTAKING**

I/we hereby certify that

- (i) All the information provided in this report are correct to the best of my knowledge and belief.
- (ii) In case, it is found, at any stage, that the information provided by me/us is false and/ or not authentic, then I/we hereby undertake to comply with the regulatory action(s) enforced against me/us and our institution, in accordance with the applicable Rules.

Place:

Signature:

Date:

Name of the Service Engineer:

Name of Supplier/Manufacturer:

Seal of Supplier/ Manufacturer:

Date:

Signature:

Name of Nuclear Medicine Technologist:

Date:

Signature:

(Licensee, Nuclear Medicine Facility)

Seal of the institute

**Document to be attached:**

- A. Manufacturer's specification sheet indicating the reference values mentioned
- B. QA report of the CT component as per the AERB format if applicable

**Reference documents:**

1. National Electrical Manufacturers Association. "NEMA Standards Publication NU 2-2001: Performance Measurements of Positron Emission Tomographs," USA, 2001.
2. National Electrical Manufacturers Association. "NEMA Standards Publication NU 2-2007: Performance Measurements of Positron Emission Tomographs," USA, 2007.
3. International Atomic Energy Agency. "IAEA Health Human Series No. 1: Quality Assurance for PET and PET/CT Systems," Vienna, Austria, 2009.

## II. Test related to Computed Tomography (CT) component

### QUALITY ASSURANCE TEST REPORT FOR COMPUTED TOMOGRAPHY EQUIPMENT

#### A. DETAILS OF THE DIAGNOSTIC X-RAY EQUIPMENT

1	Name of the Institution and City	
2	Type of Equipment	
3	Model Name	
4	Name of the Manufacturer	
5	Name(s) of Person(s) testing the equipment and Name of Supplier/Service Agency	
6	Date and Duration of the Tests	

#### B. SUMMARY OF MECHANICAL SAFETY PERFORMANCE TESTS REPORT

Sr. No.	Parameters Tested	Specific Value	Measured Value	Tolerance	Remarks
1.	Alignment of Table to Gantry			$\pm 5$ mm	
2.	Accuracy of Gantry Tilt			$\pm 2^0$	
3.	Table Indexing Accuracy			$\pm 1.0$ mm	

#### C. SUMMARY OF RADIATION SAFETY PERFORMANCE TEST REPORT

Sr. No.	Parameter Tested	Specific Value	Measured Value	Tolerance	Remarks								
1.	Slice Thickness (mm)			<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 60%;">For Slice Thickness</td> <td style="width: 40%;">0.5 mm</td> </tr> <tr> <td>a. Less than 1 mm</td> <td></td> </tr> <tr> <td>b. 1 mm to 2 mm</td> <td><math>\pm 50\%</math></td> </tr> <tr> <td>c. Above 2 mm</td> <td><math>\pm 1</math> mm</td> </tr> </table>	For Slice Thickness	0.5 mm	a. Less than 1 mm		b. 1 mm to 2 mm	$\pm 50\%$	c. Above 2 mm	$\pm 1$ mm	
For Slice Thickness	0.5 mm												
a. Less than 1 mm													
b. 1 mm to 2 mm	$\pm 50\%$												
c. Above 2 mm	$\pm 1$ mm												
2.	Accuracy of Operating Potential (kV)			$\pm 2$ kV									
3.	Total Filtration			1.5 mm Al for $kV \leq 70$ 2.0 mm Al for $kV \leq 100$ 2.5 mm Al for $kV > 100$									

4.	Accuracy of Timer			Percentage Error < 10 %	
5.	Linearity of Radiation Output (mA/ mAs Linearity)			CoL < 0.1	
6.	Reproducibility of Radiation Output			CoV < 0.05	
7.	Radiation Dose Test [Weighted Computed Tomography Dose Index (CTDI <sub>w</sub> ) for Head and Body Phantom]			± 20 % of Stated Values	
8.	Low Contrast Resolution			As per Technical Specifications Or 5.0 mm at 1% contrast	
9.	High Contrast Resolution			As per Technical Specifications Or 3.12 lp/cm	
10.	Radiation Leakage Level from X-ray Tube Housing (Measurement at maximum kVp and corresponding mA)			< 1mGy in one hour	

I hereby undertake that all the information provided above is correct and in accordance with the detailed Quality Assurance Report enclosed herewith.

Place:

Date:

Signature:

Name of the Service Engineer:

Name of Supplier/QA Agency:

Seal of Supplier/QA Agency:

#Signature of Institution's Representative:

Name of Institution:

Seal of the Institution:

*# Quality Assurance Tests Report shall be signed by Institution's Representative and duly stamped by the User's Institution.*



## Part-I MECHANICAL SAFETY TESTS FOR COMPUTED TOMOGRAPHY EQUIPMENT

### A.1 Alignment of Table to Gantry

Result :

Tolerance :  $\pm 5$  mm

### A.2. Accuracy of Gantry Tilt

Exposure Parameters : kVp: mAs:

Actual Gantry Tilt :

Measured Gantry Tilt :

Result :

Tolerance :  $\pm 2^0$

### A.3. Table Indexing Accuracy

Initial table position :

Load on couch :

Exposure parameters : kVp: mAs: Slice thickness:

Applied table increments :

Table position from reference position	1 cm	2 cm	3 cm	4 cm	5 cm
Measured					

Tolerance :  $\pm 1.0$  mm

## PART-II RADIATION SAFETY TESTS OF COMPUTED TOMOGRAPHY EQUIPMENT

### 1. Slice Thickness / Radiation Profile Width:

Exposure Parameters: kVp: mAs:

Applied Slice Thickness (mm)	Measured Density Profile Width (FWHM)	Tolerance	
		For slice thickness	
		a. Less than 1 mm	0.5 mm
		b. 1 mm to 2 mm	± 50%
		c. Above 2 mm	±1 mm

**2. Accuracy of Operating Potential:**

Set kV	mA station I	mA station II	mA station III	Average kVp

Tolerance : ± 2 kVp

**3. Total Filtration**

Total Filtration (measurement at maximum kVp):

Tolerance: 1.5 mm Al for kV ≤ 70, 2.0 mm Al for kV ≤ 100, 2.5 mm Al for kV > 100

**4. Accuracy of Timer :**

Set Time	Observed Time	Percentage Error

Tolerance: Percent Error =  $|\text{Observed value} - \text{Exact Value}| / \text{Exact value} * 100 < 10 \%$

**5. Linearity of Radiation Output (Measurement of mA / mAs linearity)**

Operating Parameters: kVp:                      Slice Thickness:

Radiation Output (μGy or mR)				μGy/mAs or mR/mAs (X)
mA / mAs	I	II	III	

$$\text{Coefficient of Linearity (CoL)} = \frac{X_{\max} - X_{\min}}{X_{\max} + X_{\min}}$$

Tolerance:  $\text{CoL} < 0.1$

## 6. Reproducibility of Radiation Output

Operating Parameters : mAs:

Slice Thickness:

Operating Potential (kVp)	Radiation Output ( $\mu\text{Gy}$ or mR)					Mean (X)	CoV
	1	2	3	4	5		

$$\text{Coefficient of Variation (CoV)} = X^{-1} [\sum (X_i - X)^2 / n - 1]^{1/2}$$

Tolerance:  $\text{CoV} < 0.05$

## 7. Radiation Dose Test

### Measurement of Weighted Computed Tomography Dose Index (CTDI<sub>w</sub>)

Use pencil ionization chamber connected to a suitable electrometer, in conjunction with a head/body phantom. Measure the dose in the axial and peripheral cavities of the phantom for the techniques specified by the manufacturer.

Operating Parameters: kVp:

mAs:

Slice Thickness:

Result:

Head

Body

Axial dose : ----- mGy/mAs                      -----mGy/mAs

Peripheral dose : ----- mGy/mAs                      -----mGy/mAs

: ----- mGy/mAs                      -----mGy/mAs

: -----mGy/mAs                      -----mGy/mAs

: -----m Gy/mAs                      -----mGy/mAs

Peripheral dose (Mean): -----mGy/mAs                      -----mGy/mAs

CTDI<sub>c</sub> : ----- mGy/mAs                      -----mGy/mAs

CTDI<sub>p (mean)</sub> : ----- mGy/mAs                      -----mGy/mAs

$$\text{Weighted CTDI (CTDI}_w) = 1/3 \text{CTDI}_c + 2/3 \text{CTDI}_p$$



Workload = 500 mA-min in one hour for measurement of tube housing leakage

$$\text{Max leakage} = \frac{500 \text{ mA-min in one hour} \times \text{Max radiation leakage level (mR/hr)}}{60 \times \text{mA used for measurement}}$$

Maximum radiation leakage from tube at 1m = ----- mR in one hour

Result: Maximum radiation leakage at 1 meter from the focus of CT Tube is ..... mGy in one hour.

**Recommended upper limit:** Radiation leakage at 1 meter from the focus of CT tube should not exceed 1mGy in one hour (115 mR in one hour).

### **11. Details of Radiation Protection Survey of the Installation**

Date of radiation protection survey:

Whether radiation survey meter used for the survey has valid calibration certificate: Yes/No

Phantom: CTDI Body Phantom

Workload of the CT facility: ----- mA-min/week

Exposure Setting:-

Applied Voltage (kV):

Tube Current (mA):

Exposure Time(s):

Location	Measured radiation level (mR/hr)
Control Console (Operator Position)	
Outside Patient Entrance Door	
Behind Windows (if applicable)	
Patient Waiting Area	

$$\text{Maximum Radiation level/week (mR/wk.)} = \frac{\text{----- mA-min/week} \times \text{Max radiation level (mR/hr)}}{60 \times \text{mA used for measurement}}$$

#### **Permissible limit**

For location of Radiation Worker: 20 mSv in a year (40 mR/week)

For Location of Member of Public: 1 mSv in a year (2 mR/week)